

KCNMB1 Rabbit pAb

货号: B15341

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	14kDa/21kDa
实测分子量	18kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,THP-1,MCF7
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 30-102 of human KCNM B1 (NP_004128.1).
序列	TYIILVTTVLPLYQKSVWTQESKCHLIETNIRDQEELKGKKVPQYPCLWVNVSAAGRWAFLYHTEDTRDQNQQ

靶点信息

研究背景	MaxiK channels are large conductance, voltage and calcium-sensitive potassium channels which are fundamental to the control of smooth muscle tone and neuronal excitability. MaxiK channels can be formed by 2 subunits: the pore-forming alpha subunit and the product of this gene, the modulatory beta subunit. Intracellular calcium regulates the physical association between the alpha and beta subunits.
基因ID	3779

基因名	KCNMB1
Swiss	Q16558
别名	KCNMB1;BKbeta1;K(VCA)beta;SLO-BETA;hbeta1;hslo-beta;k(VCA)beta-1;slo-beta-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn