

MAN2A1 Rabbit pAb

货号: B15344

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	131kDa
实测分子量	131-170kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,MCF7,NCI-H460,HT-29,HT-1080,Mouse liver,Rat liver
细胞定位	Golgi apparatus membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 865-1144 of human MA N2A1 (NP_002363.2).
序列	SNIVDIRKVYNREIAMKISSDIKSQNRFYTDLNGYQIQPRMTLSKLPLQANVYPMTTMAYIQDAKHRRLTLLSAQSLGVSSLNS GQIEVIMDRRLMQDDNRGLEQGIQDNKITANLFRILLEKRSVANTEEEEKSVSYPSLLSHITSSLMNHVPVPMANKFSSPTLE LQGEFSPLQSSLP CDIHLVNLRTIQSKVGNGHSNEAALILHRKGFDCRFSSKGTGLFCSTTQGKILVQKLLNKFIVESLTPSS LSLMHSPPGTQNISEINLSPMEISTFRIQLR

靶点信息

研究背景	This gene encodes a glycosyl hydrolase that localizes to the Golgi and catalyzes the final hydrolytic step i n the asparagine-linked oligosaccharide (N-glycan) maturation pathway. Mutations in the mouse homolog of this gene have been shown to cause a systemic autoimmune disease similar to human systemic lupus erythematosus.
------	---

基因ID	4124
基因名	MAN2A1
Swiss	Q16706
别名	MAN2A1;AMan II;GOLIM7;MANA2;MANII;AManII

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn