

# Phospho-eIF4EBP1-T45 Rabbit pAb

货号: **B15371**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:100
理论分子量	12kDa
实测分子量	18kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	cytoplasm,cytosol
纯化	Affinity purification

## 抗原信息

抗原信息	A phospho specific peptide corresponding to residues surrounding T45 of human Phospho-eIF4EBP1-T45
序列	

## 靶点信息

研究背景	This gene encodes one member of a family of translation repressor proteins. The protein directly interact s with eukaryotic translation initiation factor 4E (eIF4E), which is a limiting component of the multisubunit complex that recruits 40S ribosomal subunits to the 5' end of mRNAs. Interaction of this protein with eIF4 E inhibits complex assembly and represses translation. This protein is phosphorylated in response to vari ous signals including UV irradiation and insulin signaling, resulting in its dissociation from eIF4E and activ ation of mRNA translation.
基因ID	1978

基因名	EIF4EBP1
Swiss	Q13541
别名	4E-BP1;4EBP1;BP-1;PHAS-I;EIF4EBP1

产品验证

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)