

RBM39 Rabbit pAb

货号: B15379

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	57kDa/58kDa/59kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,U-251MG,293T,LO2,Mouse brain,Mouse heart
细胞定位	Nucleus speckle
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 400 to the C-terminus of human RBM 39 (NP_001229528.1).
序列	ATQCFQLSNMFNPQTETEEVGWDTEIKDDVIEECNKHGGVIHIYVDKNSAQGNVYVKCPSIAAAIAAVNALHGRWFAGKMI TAAYVPLPTYHNLFPDSMTATQLLVPSRR

靶点信息

研究背景	This gene encodes a member of the U2AF65 family of proteins. The encoded protein is found in the nucle us, where it co-localizes with core spliceosomal proteins. It has been shown to play a role in both steroid hormone receptor-mediated transcription and alternative splicing, and it is also a transcriptional coregula tor of the viral oncoprotein v-Rel. Multiple transcript variants have been observed for this gene. A related pseudogene has been identified on chromosome X.
------	---

基因ID	9584
基因名	RBM39
Swiss	Q14498
别名	RBM39;CAPER;CAPERalpha;FSAP59;HCC1;RNPC2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn