

MCM3AP Rabbit pAb

货号: B15383

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	80kDa/218kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa,COS-7,NIH/3T3,Mouse brain
细胞定位	Cytoplasm,Nucleus,Nucleus envelope,nuclear pore complex,nucleoplasm
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 250-550 of human MCM 3AP (NP_003897.2).
序列	SSFPVSSAVLGEPFQASKAGVRQGCCEEAVSQVEPLPSLMKGLKRKEDQDRSPRRHGHEPAEDSDPLSRGDHPPDKRPV RLNRPRGGTLFGRTIQDVFKSNKEVGRLGNKEAKKETGFVESAESDHMAIPGGNQSVLAPSRIPGVNKEEETESREKKED SLRGTPARQSNRSESTDSLGGLSPSEVTAIQCKNIPDYLNDRTILENHFGKIAKVQRIFTRRSKKLAVVHFFDHASALARK KGKSLHKDMAIFWHRKKISPNNKPFSLKEKKPGDGEVSPSTEDAPFQHSPLGKAAGRTGAS

靶点信息

研究背景	The minichromosome maintenance protein 3 (MCM3) is one of the MCM proteins essential for the initiation of DNA replication. The protein encoded by this gene is a MCM3 binding protein. It was reported to have phosphorylation-dependent DNA-primase activity, which was up-regulated in antigen immunization induced germinal center. This protein was demonstrated to be an acetyltransferase that acetylates MCM3 and plays a role in DNA replication. The mutagenesis of a nuclear localization signal of MCM3 affects the binding of this protein with MCM3, suggesting that this protein may also facilitate MCM3 nuclear localization. This gene is expressed in the brain or in neuronal tissue. An allelic variant encoding amino acid Lys at 915, instead of conserved Glu, has been identified in patients with mild intellectual disability.
基因ID	8888
基因名	MCM3AP
Swiss	O60318
别名	MCM3AP;GANP;MAP80;SAC3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn