

PRKACB Rabbit pAb

货号: B15385

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	29kDa/36kDa/39kDa/40kDa/41kDa/46kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Rat brain
细胞定位	Cell membrane,Cytoplasm,Lipid-anchor,Membrane,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-351 of human PRKACB (NP_002722.1).
序列	MGNAATAKKGSEVESVKEFLAKAKEDFLKKWENPTQNNAGLEDFERKKTLTGTSFGRVMLVKHKATEQYYAMKILDKQK VVKLKQIEHTLNEKRILQAVNFPFLVRLEYAFKDNSNLYMVM EYVPGGEMFSLRRIGRFSEPHARFYAAQIVLTFEYLHSLD LIYRDLKPENLLIDHQGYIQVTDGFAKRVKGRTWTL CGTP EYLAPEIILSKGYNKAVDWWALGVLIYEMAAGYPPFFADQPI QIYEKIVSGKVRFP SHFSSDLKDLLRNLLQVDLTKRFGNLKNGVSDIKTHKWFATTDWIAIYQRKVEAPFIPKFRGSGDTSN FDDYEEEDIRVSITEKCAKEFG EF

靶点信息

研究背景	The protein encoded by this gene is a member of the serine/threonine protein kinase family. The encoded protein is a catalytic subunit of cAMP (cyclic AMP)-dependent protein kinase, which mediates signalling through cAMP. cAMP signaling is important to a number of processes, including cell proliferation and differentiation. Multiple alternatively spliced transcript variants encoding distinct isoforms have been observed.
------	---

基因ID	5567
基因名	PRKACB
Swiss	P22694
别名	PRKACB;PKA C-beta;PKACB;PKAC-beta

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn