

ESYT1 Rabbit pAb

货号: B15414

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:200 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	122kDa/124kDa
实测分子量	123kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,293T,LO2,Jurkat,HeLa,Mouse brain,Mouse heart
细胞定位	Cell membrane,Endoplasmic reticulum membrane,Multi-pass membrane protein,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 90-220 of human ESYT1 (NP_056107.1).
序列	GWRRVRDEKERSLRAARQLLDDEEQLTAKTLYMSHRELPWVSFPDVEKAEWLNKIVAQVWPFLGQYMEKLLAETVAPAVRGSNPHLQFTFTRVELGEKPLRIIGVKVHPGQRKEQILLDLNISYVGDVQ

靶点信息

研究背景	Binds glycerophospholipids in a barrel-like domain and may play a role in cellular lipid transport (By similarity). Binds calcium (via the C2 domains and translocates to sites of contact between the endoplasmic reticulum and the cell membrane in response to increased cytosolic calcium levels. Helps tether the endoplasmic reticulum to the cell membrane and promotes the formation of appositions between the endoplasmic reticulum and the cell membrane.
------	--

基因ID	23344
基因名	ESYT1
Swiss	Q9BSJ8
别名	ESYT1;FAM62A;MBC2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn