

THEMIS2 Rabbit pAb

货号: B15421

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:200 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	13kDa/29kDa/50kDa/56kDa/57kDa/72kDa
实测分子量	72kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,MCF7,Mouse spleen
细胞定位	cytoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human THEMIS 2 (NP_004839.2).
序列	MEPVPLQDFVRALDPASLPRVLRVCSGVYFEGSIYEISGNECCLSTGDLIKVTQVRLQKVVCEPNKTSQTMELAPNFQGYFT PLNTPQSYETLEELVSATTQSSKQLPTCFMSTHRIVTEGRVVTEDQLLMLEAVVMHLGIRSARCVLGMEGQQVILHPLSQK GPFWTWEPAPRTLLQVLQDPALKDLVLTCTLPWHSILIRPQYEIQAIMHIFSSLRIAATRSAQAQTQGEDLARVHQGWLQ YVQQDSCPQEGPQAR

靶点信息

研究背景	May constitute a control point in macrophage inflammatory response, promoting LPS-induced TLR4-mediated TNF production. Determines the threshold for activation of B cells by low-affinity and low-avidity ligands via PLCG2 activation and its downstream pathways (By similarity).
------	--

基因ID	9473
基因名	THEMIS2
Swiss	Q5TEJ8
别名	THEMIS2;C1orf38;ICB-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn