

GPR87 Rabbit pAb

货号: B15432

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	41kDa
实测分子量	41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	DU145
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human GPR87 (NP_076404.3).
序列	HNQSIRVVAVFFTCFLPYHLCRIPFTFSLDRLLDESAQKILYYCKEITLFLSACNVCLDPIIYFFMCRSFSRRLFKKSNIRTRS ESIRSLQSVRRSEVR

靶点信息

研究背景	This gene encodes a G protein-coupled receptor and is located in a cluster of G protein-couple receptor genes on chromosome 3. The encoded protein has been shown to be overexpressed in lung squamous cell carcinoma (PMID:18057535) and regulated by p53 (PMID:19602589).
基因ID	53836

基因名	GPR87
Swiss	Q9BY21
别名	GPR87;FKSG78;GPR95;KPG_002

产品验证



Western blot analysis of GPR87 expressed in DU145 using GPR87 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn