

# VPS36 Rabbit pAb

货号: B15433

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IF/ICC                                                                                            |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IF/ICC:</b> 1:50 - 1:200                                             |
| 理论分子量 | 36kDa/43kDa                                                                                          |
| 实测分子量 | 43kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Rat liver,Rat kidney                                                                                 |
| 细胞定位  | Cytoplasm,Endosome,Late endosome,Membrane,Nucleus                                                    |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                           |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 70-220 of human VPS36 (NP_057159.2).                                        |
| 序列   | EEQAAGIGKSAKIVVHLHPAPPNKEPGPFQSSKNSYIKLSFKEHGQIEFYRRLSEEMTQRRWENMPVSQSLQTNRGPQPG RIRAVGIVGIERKLEEKRKETDKNISEAFEDLSKLMIAKAKEMVELSKSIANKIKDKQGDITEDETIRFKSY |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                 |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a protein that is a subunit of the endosomal sorting complex required for transport II ( ESCRT-II). This protein complex functions in sorting of ubiquitinated membrane proteins during endocytosis. A similar protein complex in rat is associated with RNA polymerase elongation factor II. |
| 基因ID | 51028                                                                                                                                                                                                                                                                                                           |

|       |                             |
|-------|-----------------------------|
| 基因名   | VPS36                       |
| Swiss | Q86VN1                      |
| 别名    | VPS36;C13orf9;CGI-145;EAP45 |

产品验证



Western blot analysis of VPS36 expressed in Rat liver,Rat kidney using VPS36 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)