

RASAL2 Rabbit pAb

货号: B15435

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	128kDa/143kDa
实测分子量	160kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 750-1040 of human RAS AL2 (NP_733793.2).
序列	LADITKSLTNPTPIQQQLRRFTEHNSSPNVSGSLSSGLQKIFEDPTDSLHKLKSPSQDNTDSYFRGKTLVVQQASSQSM TYSEKDERESSLPNGRSVSLMDLQDTHAAQVEHASVMLDVPIRLTGSQLSITQVASIKQLRETQSTPQSAPQVRRPLHPAL NQPGGLQPLSFQNPVYHLNNPAMPKASIDSSLENLSTASSRSQSNSEDFKLSGSPSNSSMEDFTKRSTQSEDFSRRHV PDRHIPLALPRQNSTGQAQIRKVDQGGLGARAKAPPSLPHSASLRSTGS

靶点信息

研究背景	This gene encodes a protein that contains the GAP-related domain (GRD), a characteristic domain of GTP ase-activating proteins (GAPs). GAPs function as activators of Ras superfamily of small GTPases. The prot ein encoded by this gene is able to complement the defective RasGAP function in a yeast system. Two alt ernatively spliced transcript variants of this gene encoding distinct isoforms have been reported.
------	---

基因ID	9462
基因名	RASAL2
Swiss	Q9UJF2
别名	RASAL2;NGAP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn