

HOXC13 Rabbit pAb

货号: B15443

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	35kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human HOXC13 (NP_059106.2).
序列	MTTSLLLHPRWPESLMYVYEDSAAESGIGGGGGGGGGGTGGAGGGCSGASPGKAPSM DGLGSSCPASHCRDLLPHPV LGRPPAPLGAPQGAVYTDIPAPEAARQCAPPAPPTSSSATLGYGYPFGGSYYGCRLSHNVNLQKPCAYHPGDKYPEPS GALPGDDLSSRAKEFAFYPSFASSYQAMPGYLDVSVVPGISGHPEPRHDALIPVEGYQHWALSNGWDSQVYCSKEQSQS AHLWKSPFPDVVPLQPEVSSYRRGRKKRVPTYKVQLKELEKEYAASKFITKEKRRRISATTNLSERQVTIWFQNRrvKEKKV VSKSKAPHLHST

靶点信息

研究背景	This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, which are located on different chromosomes and consist of 9 to 11 genes arranged in tandem. This gene is one of several homeobox HOXC genes located in a cluster on chromosome 12. The product of this gene may play a role in the development of hair, nail, and filiform papilla.
基因ID	3229
基因名	HOXC13
Swiss	P31276
别名	HOXC13;ECTD9;HOX3;HOX3G

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn