

ACTL8 Rabbit pAb

货号: B15450

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	41kDa
实测分子量	41-50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2,Mouse testis,Mouse pancreas,Mouse brain
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-366 of human ACTL8 (NP_110439.2).
序列	MAARTVIIDHGSGFLKAGTAGWNEPQMVFNPVNYLPCKENPGPSYARRRVSLGIDICHPTFSYPIERGRILNWEGVQYL WSFVLENHRREQEVPPVIITETPLREPADRKKMLEILFELLHVPSVLLADQLQMSLYASGLLTGVVVDSDGYGLTRVQPFHQG RPLPASGKTLEFAGQDLSAYLLKSLFKEDCDRRCLFQLETVAVTQMKNKYVPQNLGEALDFRERQQSALDESNTYQLPDG SRVELTPMQRVAPEMFFSPQVFEQGPSIPRAIVESVESCEISLRPLLVSHVMACGGNTLYPGFTKRLFRELMGDHVSSTKA TVWEGSNRNFVWLGASVVAHLSTYQSEWMSREEYGEHMRM

靶点信息

研究背景	Involved in epithelial cell differentiation. Predicted to be located in cytoplasm. Predicted to be part of dyn actin complex.
基因ID	81569

基因名	ACTL8
Swiss	Q9H568
别名	ACTL8;CT57

产品验证



Western blot analysis of ACTL8 expressed in Mouse testis,HeLa,HepG2 using ACTL8 Rabbit pAb at 1:100 0. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Bloc king buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn