

ZCCHC17 Rabbit pAb

货号: B15454

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	24kDa/26kDa/27kDa
实测分子量	27kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-241 of human ZCCHC 17 (NP_057589.2).
序列	MNSGRPETMENLPALYTIFQGEVAMVTDYGAFIKIPGCRKQGLVHRTHMSSCRVDKPSEIVDVGDKVWVKLIGREMKNDR IKVSLSMKVVNQGTGKDLPNNVIIIEQEERRRSFQDYTGQKITLEAVLNTTCKKCGCKGHFAKDCFMQPGGTKYSLIPDE EEEKEEAKSAEFEKDPTRNPSRKRKKEKKKKKHRDRKSSDSDSSDSES DTGKRARHTSKDSKAAKKKKKKKKHKKKKHK E

靶点信息

研究背景	Enables identical protein binding activity. Predicted to be involved in RNA stabilization. Predicted to be lo cated in nucleolus. Predicted to be part of cytosolic large ribosomal subunit.
基因ID	51538

基因名	ZCCHC17
Swiss	Q9NP64
别名	ZCCHC17;HSPC251;PS1D;pNO40

产品验证



Western blot analysis of ZCCHC17 expressed in Mouse brain using ZCCHC17 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn