

ELP5 Rabbit pAb

货号: B15456

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	19kDa/30kDa/34kDa/35kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,LO2,A375,Mouse brain,Mouse skeletal muscle,Mouse testis
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 17-316 of human ELP5 (NP_056177.3).
序列	MLDSLALGGLVLLRDSVEWEGRSLLKALVKKSALCGEQVHILGCEVSEEEFREGFDSDINNRLVYHDFFRDPLNWSKTEE AFPGGPLGALRAMCKRTDPVPVTIALDSLWLLLRPCTTLCQVLHAVSHQDSCPGDSSSVGKVSVLGLLHEELHGPGPV GALSSLAQTEVTLGGMGQASAHILCRRPRQRPTDQTQWFSILPDFSLDLQEGPSVESQPYSDPHIPPPVDPTTHLTFNLHL SKKEREARDSLILPFQFSSEKQKALLRPRPGQATSHIFYEPDAYDDLDQEDPDDDLDI

靶点信息

研究背景	Component of the RNA polymerase II elongator complex, a multiprotein complex associated with the RNA polymerase II (Pol II holoenzyme, and which is involved in transcriptional elongation. The elongator complex catalyzes formation of carboxymethyluridine in the wobble base at position 34 in tRNAs. Involved in cell migration (By similarity. May be involved in TP53-mediated transcriptional regulation.
------	---

基因ID	23587
基因名	ELP5
Swiss	Q8TE02
别名	ELP5;C17orf81;DERP6;HSPC002;MST071;MSTP071

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn