

JRKL Rabbit pAb

货号: B15463

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	59kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,SH-SY5Y,Mouse brain,Mouse thymus
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 225-524 of human JRKL (NP_003763.2).
序列	LCVVGKAKKPRSFKSTDTLNLPSYFSQKGAWMDLSIFRQWFDKIFVPQVREYLRSKGLQEKAVALLLDNSPTHPNENVLR SDDGQIFAKYLPPNVASLIQPSDQGVIATMKRNYRAGLLQNNLEEGNDLKSFWKKLTLLDALYEIAMAWNLVKPTISRAW KKILPMVEEKESLDFDVEDISVATVAAILQHTKGLENVTTENLEKWLEV DSTPGYEVLT DSEIIRRAQQGADESSNEEEEEIE L IPEKHINHAAALQWTENLLDYLEQQGDMILPDRLVIRKLRATIRNKQKMTKSSQ

靶点信息

研究背景	The function of this gene has not yet been defined, however, the encoded protein shares similarity with t he human (41% identical) and mouse (34% identical) jerky gene products. This protein may act as a nucl ear regulatory protein. Alternatively spliced transcript variants have been found for this gene.
基因ID	8690

基因名	JRKL
Swiss	Q9Y4A0
别名	JRKL;HHMJG;JRK-like

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn