

ZNF22 Rabbit pAb

货号: B15465

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	25kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-224 of human ZNF22 (NP_008894.2).
序列	MRLAKPKAGISRSSSQGKAYENKRKTGRQRQKWGMTIRFDSSFSRLRRSLDDKPYKCTECEKSFSQSSTLFQHQKIHTGK KSHKCADCGKSFFQSSNLIQHRRihtGEKPYKCDECGESFKQSSNLIQHqRIHTGEKPYQCDECGRCFSQSSHLIQHQRT HTGEKPYQCSECGKCFsQSSHLRQHMKVHKEEKPRKTRGKNIRVKThLPsWKAGTGRKSVAGLR

靶点信息

研究背景	Binds DNA through the consensus sequence 5'-CAATG-3'. May be involved in transcriptional regulation and may play a role in tooth formation (By similarity).
基因ID	7570

基因名	ZNF22
Swiss	P17026
别名	ZNF22;HKR-T1;KOX15;ZNF422;Zfp422

产品验证



Western blot analysis of ZNF22 expressed in Mouse liver using ZNF22 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn