

GPR6 Rabbit pAb

货号: B15468

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	37kDa/39kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 253-362 of human GPR6 (NP_005275.1).
序列	WRHAHQIALQQHCLAPPHLAATRKGVGTLAVVLGTFGASWLPFAIYCVVGSHPEDPAVYTYATLLPATYNSMINPIIYAFRNQ EIQRALWLLLCGCFQSKVPFRSRSPSEV

靶点信息

研究背景	Orphan receptor with constitutive G(s signaling activity that activate cyclic AMP. Promotes neurite outgro wth and blocks myelin inhibition in neurons (By similarity.
基因ID	2830

基因名	GPR6
Swiss	P46095
别名	GPR6

产品验证



Western blot analysis of GPR6 expressed in U-251MG using GPR6 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn