

ETNK2 Rabbit pAb

货号: B15471

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	40kDa/44kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,DU145,HepG2
细胞定位	cytoplasm,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-386 of human ETNK2 (NP_060678.2).
序列	MAVPPSAPQPRASFHLRRHTPCPQCSWGMEEKAAASASCREPPGPPRAAAVAYFGISVDPDDILPGALRLIQELRPHWKPE QVRTKRFTDGITNKLVCYVEEDMQDCVLVRVYGERTELLVDRENEVRNFQLLRAHSCAPKLYCTFQNGLCYEYMQGVAL EPEHIREPRLFRLIALEMAKIHTIHANGSLPKPILWHKMHNYYFTLVKNEINPSLSADVPKVEVLERELAWLKEHLSQLESPVVF CHNDLLCKNIIYDSIKGHVRFIDYEYAGYNYQAFDIGNHFNFAFVNEVDYCLYPARETQLQWLHYYLQAQKGMVTPREV QRLYVQVNFALASHFFWALWALIQNQYSTIDFDFLRYAVIRFNQYFKVKPQASALEMPK

靶点信息

研究背景	The protein encoded by this gene is a member of choline/ethanolamine kinase family which catalyzes the first step of phosphatidylethanolamine (PtdEtn) biosynthesis via the cytidine diphosphate (CDP) ethanolamine pathway. Alternative splicing results in multiple transcript variants.
------	--

基因ID	55224
基因名	ETNK2
Swiss	Q9NVF9
别名	ETNK2;EKI2;HMFT1716

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn