

TMTC3 Rabbit pAb

货号: B15476

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	103kDa/104kDa
实测分子量	104kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HT-29
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 615-914 of human TMT C3 (NP_861448.2).
序列	EALKNFNRALELNPKHKLALFNSAIVMQESGEVKLRPEARKRLLSYINEEPLDANGYFNLGMLAMDDKKDNEAEIWMKKAI KLQADFRSALFNLALLYSQTAKELKALPILELLRYPDHIKGLILKGDILMNQKKDILGAKKCFERILEMDPSNVQGKHNL CV VYFEEKDLLKAERCCLLETALAPHEEYIQRHLNIVRDKISSSFIEPIFPTSKISSVEGKKIPTESVKEIRGESRQTQIVKTS DNKS QSKSNKQLGKNGDEETPHKTTKDIKEIEKKRVAALKRLEEIERILNGE

靶点信息

研究背景	This gene encodes a protein that belongs to the transmembrane and tetratricopeptide repeat-containing protein family.
基因ID	160418

基因名	TMTC3
Swiss	Q6ZXV5
别名	TMTC3;LIS8;SMILE

产品验证



Western blot analysis of TMTC3 expressed in HeLa,HT-29 using TMTC3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5 % non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn