

# FSHB Rabbit pAb

货号: B15479

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Mouse,Rat                                                                                            |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 14kDa                                                                                                |
| 实测分子量 | 14kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Mouse brain,Mouse lung,Rat brain                                                                     |
| 细胞定位  | Secreted                                                                                             |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 19-129 of human FSHB ( NP_000501.1). |
| 序列   | NSCELTNITIAIEKEECRFCISINTTWCAGYCYTRDLVYKDPARPKIQKTCTFKELVYETVRVPGCAHHADSLYTPVATQCHC<br>GKCDSDDTDCTVRGLGPSYCSFGEMKE |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The pituitary glycoprotein hormone family includes follicle-stimulating hormone, luteinizing hormone, cho<br>rionic gonadotropin, and thyroid-stimulating hormone. All of these glycoproteins consist of an identical alp<br>ha subunit and a hormone-specific beta subunit. This gene encodes the beta subunit of follicle-stimulatin<br>g hormone. In conjunction with luteinizing hormone, follicle-stimulating hormone induces egg and sperm<br>production. Alternative splicing results in two transcript variants encoding the same protein. |
| 基因ID | 2488                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |

|       |           |
|-------|-----------|
| 基因名   | FSHB      |
| Swiss | P01225    |
| 别名    | FSHB;HH24 |

产品验证



Western blot analysis of FSHB expressed in Mouse brain,Mouse lung,Rat brain using FSHB Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)