

CCDC83 Rabbit pAb

货号: B15486

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	37kDa/48kDa/52kDa
实测分子量	49kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,HeLa
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CCDC83 (NP_001273088.1).
序列	MENSGKANKKDTHDGPPKEIKLPTSEALLDYQCQIKEDAVEQFMFQIKTLRKKNQKYHERNSRLKEEQIWHIRHLLKELSEE KAEGLPVVTREDVEEAMKEKWKFERDQEKNLDRMRMQISNAEKLFLLEKLSEKEYWEEYKNVGSERHAKLITSLQNDINTV KENAEKMSEHYKITLEDTRKKIIKETLLQLDQKKEWATQNAVKLIDKGSYLEIWENDWLKKEVAIHRKEVEELKNAIHELEAE NLVLIDQLSNCRLVDLKIPRRLYLTAAGLEVPPEEMSLELPETHIEEKSELQPT

靶点信息

研究背景	CCDC83 (Coiled-Coil Domain Containing 83) is a Protein Coding gene. Diseases associated with CCDC83 include Acute Hemorrhagic Conjunctivitis.
基因ID	220047

基因名	CCDC83
Swiss	Q8IWF9
别名	CCDC83;CT148;HSD9

产品验证



Western blot analysis of CCDC83 expressed in A-431,HeLa using CCDC83 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn