

# UBQLNL Rabbit pAb

货号: B15487

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IF/ICC                                                                                            |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IF/ICC:</b> 1:50 - 1:200                                             |
| 理论分子量 | 29kDa/52kDa                                                                                          |
| 实测分子量 | 53kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Mouse heart,Rat heart                                                                                |
| 细胞定位  | cytosol                                                                                              |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                                                                                                                                             |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 176-475 of human UBQ LNL (NP_659490.4).                                                                                                                                                                                       |
| 序列   | LENPSIQRLLSNMEFMWQFISEHLDLTQQLMQQNPEVSRLLLDNSEILLQTLLEARNLAMIQEIMQIQPSQNLEYPLNPQPY LGLETMPGGNNALGQNYADINDQMLNSMQDPFGGNPFTALLAGQVLEQVQSSPPPPPSQEQDQLTQHPATRVYNS SGGFSSNTSANDTLNKVNHTSKANTAMISTKGQSHICATRQPAWIPALPSIELTQQLQEEYKDATVSLSSSRQTLKGDQLQ SDEQSSSQITGGMMQLLMNNPYLAAQIMLFTSMPLSEQWRQQLPTFLQQTQISDLLSA |

靶点信息

|      |                                                                                                                                                                                                   |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | Predicted to enable polyubiquitin modification-dependent protein binding activity. Predicted to be involve d in ubiquitin-dependent protein catabolic process. Predicted to be active in cytosol. |
| 基因ID | 143630                                                                                                                                                                                            |

|       |        |
|-------|--------|
| 基因名   | UBQLNL |
| Swiss | Q8IYU4 |
| 别名    | UBQLNL |

产品验证



Western blot analysis of UBQLNL expressed in Mouse heart,Rat heart using UBQLNL Rabbit pAb at 1:1000 . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)