

QPRTRabbit pAb

货号: B15495

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	30kDa
实测分子量	31kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Rat liver
细胞定位	cytoplasm,cytosol,extracellular exosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-297 of human QPRT (NP_055113.2).
序列	MDAEGLALLLPVTLAALVDSWLREDCPGLNYAALVSGAGPSQAALWAKSPGVLAGQPFFDAIFTQLNCQVSWFLPEGS KLVPVARVAEVRGPAHCLLLGERVALNTLARCSGIASAAAAAVEAARGAGWTGHVAGTRKTTPGFRLVEKYGLLVGGAAS HRYDLGGLVMVKDNHVVAAGGVEKAVRAARQAADFALKVEVECSSLQEAVQAAEAGADLVLLDNFKPEELHPTATVLKA QFPSVAVEASGGITLDNLPQFCGPHIDVISMGMILTQAAPALDFSLKLFakeVAPVPKIH

靶点信息

研究背景	This gene encodes a key enzyme in catabolism of quinolinate, an intermediate in the tryptophan-nicotina mide adenine dinucleotide pathway. Quinolinate acts as a most potent endogenous exitotoxin to neurons . Elevation of quinolinate levels in the brain has been linked to the pathogenesis of neurodegenerative dis orders such as epilepsy, Alzheimer's disease, and Huntington's disease. Alternative splicing results in mul tiple transcript variants.
------	--

基因ID	23475
基因名	QPRT
Swiss	Q15274
别名	QPRT;HEL-S-90n;QPRTase

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn