

ADAM23 Rabbit pAb

货号: B15500

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	91kDa/92kDa
实测分子量	92kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Raji,Mouse heart
细胞定位	Cell membrane,Secreted,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 350-500 of human ADA M23 (NP_003803.1).
序列	VETWTEKDQIDITTNPVQMLHEFSKYRQRIQHADAVHLISRVTFFHYKRSSLSYFGGVCSRTRGVGVNEYGLPMAVAQVLSQSLAQNLGIQWEPSSRKPKCDCTESWGGCIMEETGVSHSRKFSKCSILEYRDFLQRGGGACLFNRPTKLF

靶点信息

研究背景	This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. It is reported that inactivation of this gene is associated with tumorigenesis in human cancers.
基因ID	8745

基因名	ADAM23
Swiss	O75077
别名	ADAM23;MDC-3;MDC3

产品验证



Western blot analysis of ADAM23 expressed in 293T,Raji,Mouse heart using ADAM23 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn