

KIF19 Rabbit pAb

货号: B15519

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	57kDa/61kDa/111kDa
实测分子量	111kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Mouse testis,Mouse lung,Rat lung
细胞定位	Cell projection,Cytoplasm,cilium,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human KIF19 (NP_694941.2).
序列	MKDSGDSKDQQLMVALRVRPISVAELEEATLIAHKVDEQMVVLMDPMEDPDDILRAHRSREKSYLFDVAFDFTATQEM VYQATTKSLIEGVISGYNATVFAYGPTGCGKTYTMLGTDQEPGIYVQTLNDLFRAIEETSNDMEYEVMSYLEIYNEMIRDLL NPSLGYLELREDSKGVIVAGITEVSTINAKEIMQLLMKGNRQRTQEPTAANQTSSRSHAVLQVTVRQSRVKNILQEV RQ GRLFMIDLAGSERASQTQNRGQRMKEGAHINRSLALGNCINALSDKGSNKYINYRD

靶点信息

研究背景	Plus end-directed microtubule-dependent motor protein that regulates the length of motile cilia by mediating depolymerization of microtubules at ciliary tips.
基因ID	124602

基因名	KIF19
Swiss	Q2TAC6
别名	KIF19;KIF19A

产品验证



Western blot analysis of KIF19 expressed in Mouse testis,Rat lung,U-87MG using KIF19 Rabbit pAb at 1: 1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn