

GPIHBP1 Rabbit pAb

货号: B15542

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	19kDa
实测分子量	20kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,LO2
细胞定位	Apical cell membrane,Basolateral cell membrane,Cell membrane,Extracellular side,GPI-anchor,Lipid-anchor,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 50-151 of human GPIHB P1 (NP_835466.2).
序列	ETNRLPGGRSRVLLRCYTCKSLPRDERCNLTQNC SHGQTCTTLIAHGNTESGLLTHTSTWCTDSCQPITKTVEGTQVTMTC CQSSLCNVPPWQSSRVQDPTG

靶点信息

研究背景	This gene encodes a capillary endothelial cell protein that facilitates the lipolytic processing of triglycerid e-rich lipoproteins. The encoded protein is a glycosylphosphatidylinositol-anchored protein that is a mem ber of the lymphocyte antigen 6 (Ly6) family. This protein plays a major role in transporting lipoprotein lip ase (LPL) from the subendothelial spaces to the capillary lumen. Mutations in this gene are the cause of h yperlipoproteinemia, type 1D. Alternate splicing results in multiple transcript variants.
------	--

基因ID	338328
基因名	GPIHBP1
Swiss	Q8IV16
别名	GPIHBP1;GPI-HBP1;HYPL1D

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn