

SRBD1 Rabbit pAb

货号: B15543

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	69kDa/111kDa
实测分子量	80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,A-549,Mouse testis,Rat testis,Rat brain
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 330-600 of human SRB D1 (NP_060549.4).
序列	LGLEGAARALLEKPGELSLLSYIRPDVKGLSTLQDIEIGVQHILADMIAKDKDTLDFIRNLCQKRHVCIQSSLAKVSSKKVNEK DVDKFLLYQHFSNIRNIHHHQILAINRGENLKVLTVKVNISDGVKDEFCEWCIQNRWRPRSFARPELMKILYNSLNSDFKR LIYPLLCREFRAKLTSDAEKESVMMFGRNLRQLLLTSPVPGRTLGMGVDPGYKHGCKLAIISPTSQILHTDVVYLHCGQGFRE AEIKITLLNFNCSTVVVINGTA

靶点信息

研究背景	Predicted to enable mRNA binding activity. Predicted to be a structural constituent of ribosome. Predicted to be involved in translation.
基因ID	55133

基因名	SRBD1
Swiss	Q8N5C6
别名	SRBD1

产品验证



Western blot analysis of SRBD1 expressed in Mouse testis,Rat testis,293T using SRBD1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn