

# ERCC6L2 Rabbit pAb

货号: B15551

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                          |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 60kDa/81kDa/177kDa                                                                                   |
| 实测分子量 | 81kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | U-87MG,HeLa,LO2,MCF-7,Mouse testis,Mouse kidney                                                      |
| 细胞定位  | Cytoplasm,Mitochondrion,Nucleus,centrosome,cytoskeleton,microtubule organizing center                |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                                                                                                                                                          |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1262-1561 of human ERCC6L2 (NP_064592.2).                                                                                                                                                                                                  |
| 序列   | TSEERQKMLRDFYASQYPEVKEFFVDSVSQFNNSSFEGEQRTRKKSDDKRESLIKPRLSDSSETLSFKDSTNKISQVCSLKT<br>YKRKSVKFNHISYREEVFFNDAETKKSPVSSTQEIDSGKNSQASEDTVTSRSLNSESETRERRELENTMKDQQDLTRTGIS<br>RKEPLLKLENKKIENPVLENTSVISLLGDTSLDDLFKSHGNSPTQLPKKVLSGPMEKAKQRPKDFWDILNEQNDESLSKLT<br>DLAVIETLCEKAPLAAPFKRREEPATSLWKSNEKFLWKKFSPSDTDENATNTQSTT |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                               |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the Snf2 family of helicase-like proteins. The encoded protein may play a role in DNA repair and mitochondrial function. Mutations in this gene have been associated with bone marrow failure syndrome 2. Alternatively spliced transcript variants that encode different protein isoforms have been described. |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                     |
|-------|-------------------------------------|
| 基因ID  | 375748                              |
| 基因名   | ERCC6L2                             |
| Swiss | Q5T890                              |
| 别名    | ERCC6L2;BMFS2;C9orf102;RAD26L;SR278 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)