

SCGB3A2 Rabbit pAb

货号: B15553

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	10kDa
实测分子量	14kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Raji,mouse lung,rat lung
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-93 of human SCGB3A2 (NP_473364.1).
序列	MKLVTFILLVTISLCSYSATAFLINKVPLPVDKLAPLPDNLFPMDPLKLLKTLGISVEHLVEGLRKCVNELGPEASEAVKKLLEALSHLV

靶点信息

研究背景	The protein encoded by this gene is a secreted lung surfactant protein and a downstream target of thyroid transcription factor. A single nucleotide polymorphism in the promoter of this gene results in susceptibility to asthma.
基因ID	117156

基因名	SCGB3A2
Swiss	Q96PL1
别名	SCGB3A2;LU103;PNSP1;UGRP1;pnSP-1

产品验证



Western blot analysis of SCGB3A2 expressed in mouse lung, rat lung, A-549 using SCGB3A2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn