

# CD62P/P-selectin Rabbit pAb

货号: B15589

产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Mouse,Rat                                                                                              |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  |                                                                                                        |
| 应用    | WB                                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                              |
| 理论分子量 | 90kDa                                                                                                  |
| 实测分子量 | 100kDa                                                                                                 |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | Mouse kidney,Rat brain                                                                                 |
| 细胞定位  | Membrane,Single-pass type I membrane protein                                                           |
| 纯化    | Affinity purification                                                                                  |

抗原信息


|      |                                                                                                                                |
|------|--------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 220-319 of human CD62P/P-selectin (NP_002996.2). |
| 序列   | FSFNSQCSFHCTDGYQVNGPSKLECLASGIWTKPPQCLAAQCPPLKIPERGNMTCLHSAKAFQHQS SCSFSCEE GFALV GPEVVQCTASGVWTAPAPVC                         |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a 140 kDa protein that is stored in the alpha-granules of platelets and Weibel-Palade bodies of endothelial cells. This protein redistributes to the plasma membrane during platelet activation and degranulation and mediates the interaction of activated endothelial cells or platelets with leukocytes. The membrane protein is a calcium-dependent receptor that binds to sialylated forms of Lewis blood group carbohydrate antigens on neutrophils and monocytes. Alternative splice variants may occur but are not well documented. |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                                |
|-------|------------------------------------------------|
| 基因ID  | 6403                                           |
| 基因名   | SELP                                           |
| Swiss | P16109                                         |
| 别名    | SELP;CD62;CD62P;GMP140;GRMP;LECAM3;PADGEM;PSEL |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)