

GLMN Rabbit pAb

货号: B15604

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	48kDa/68kDa
实测分子量	68kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,U-87MG,HepG2,Mouse testis
细胞定位	cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 405-594 of human GLM N (NP_444504.1).
序列	RCLLNTSNHSGVEAFIIQNIKNQIDMSLKRTNNKWFTGPQLISLLDLVFLPEGAETDLLQNSDRIMASLNLLRYLVIKDNE NDNQTGLWTELGNIENNFLKPLHIGLNMSKAHYEAEIKNSQEAQKSKDLC SITVSGEEIPNMPPEMQLKVLHSALFTFDLIE SVLARVEELIEIKTKSTSEENIGIK

靶点信息

研究背景	This gene encodes a phosphorylated protein that is a member of a Skp1-Cullin-F-box-like complex. The protein is essential for normal development of the vasculature and mutations in this gene have been associated with glomuvenous malformations, also called glomangiomas. Multiple splice variants encoding different isoforms have been found for this gene.
基因ID	11146

基因名	GLMN
Swiss	Q92990
别名	GLMN;FAP;FAP48;FAP68;FKBPAP;GLML;GVM;VMGLOM;glomulin

产品验证



Western blot analysis of GLMN expressed in Mouse testis,A-549,U-87MG using GLMN Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Bl ocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn