

ADGRB3 Rabbit pAb

货号: B15607

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	55kDa/171kDa
实测分子量	180kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 70-280 of human ADGR B3 (NP_001695.1).
序列	YLKFSKKDLSCSNFSLAYQFDHFSHEKIKDLLRKNHSIMQLCNSKNAFVFLQYDKNFIQIRRVFPTNFPGLQKKGEEDQKS FFEFLVLNKVSPSQFGCHVLCTWLESCLKSENGRTESCGIMYTKCTCPQHLGEWGIDDQSLILLNNVVLPLNEQTEGCLTQ ELQTTQVCNLTREAKRPPKEEFGMMGDHTIKSQRPRSVHEKRVPEQA

靶点信息

研究背景	This p53-target gene encodes a brain-specific angiogenesis inhibitor, a seven-span transmembrane prote in, and is thought to be a member of the secretin receptor family. Brain-specific angiogenesis proteins BA I2 and BAI3 are similar to BAI1 in structure, have similar tissue specificities, and may also play a role in an giogenesis.
基因ID	577

基因名	ADGRB3
Swiss	O60242
别名	ADGRB3;BAI3

产品验证



Western blot analysis of ADGRB3 expressed in 293T using ADGRB3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn