

货号: **B15669**


反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	
实测分子量	40KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa,Mouse brain,Mouse testis
细胞定位	cytosol,nuclear speck,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-357 of human SNRNP 40 (NP_004805.2).
序列	MIEQQRKRGPELPLVPVKRQRHELLLGAGSGPGAGQQQATPGALLQAGPPRCSSLQAPIMLLSGHEGEVYCKFHPNGST LASAGFDRILLWNVYGDCDNYATLKGHSGAVMELHYNTDGSMLFSASTDKTVAVWDSETGERVKRLKGHTSFVN PARRGPQLVCTGSSDDGTVKLWDIRKKAIIQTQNTYQVLAVTFNDTSDQIISSGGIDNDIKVWDLRQNKLT YTMRGHADS VTGLSLSSEGSYLLSNAMDNTVRVWDVRPFAPKERCVKIFQGNVHNFEKNLLRCWSPDGSKIAAGSADR FVYVWDDTTS RRILYKLPGHAGSINEVAFHPDEPIIIISASSDKRLYMGEIQ

研究背景	This gene encodes a component of the U5 small nuclear ribonucleoprotein (snRNP) particle. The U5 snRNP is part of the spliceosome, a multiprotein complex that catalyzes the removal of introns from pre-messenger RNAs. [provided by RefSeq, Jul 2008]
------	---

基因ID	9410
基因名	SNRNP40
Swiss	Q96DI7
别名	40K; HPRP8BP; PRP8BP; PRPF8BP; SPF38; WDR57

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn