

FTSJ3 Rabbit pAb

货号: B15673

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	96kDa
实测分子量	125kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A375,HeLa,MCF7,HT-1080
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 608-847 of human FTSJ3 (NP_060117.3).
序列	EGEEKDGISDSDSSTSSEEEESWEPLRGKKRSRGPKSDDDGFEIVPIEDPAKHRILDPEGLALGAVIASSKKAKRDLIDNSFNRYTFNEDEGELPEWFVQEEKQHRIRQLPVGKKEVEHYRKRWREINARPIKKVAEAKARKKRRMLKRLEQTRKKAEEAVVNTVDISEREKVAQLRSLYKKAGLGKEKRHVTVVAKKGVGRKVRPAGVRGHFKVVDSDRMKKDQRAQQRKEQKKKHKRK

靶点信息

研究背景	Although the function of this gene is not known, the existence of this gene is supported by mRNA and EST data. A possible function of the encoded protein can be inferred from amino acid sequence similarity to the E.coli Ftsj protein and to a mouse protein possibly involved in embryogenesis.
基因ID	117246

基因名	FTSJ3
Swiss	Q8IY81
别名	FTSJ3;EPCS3;SPB1

产品验证



Western blot analysis of FTSJ3 expressed in A375,HeLa,MCF7 using FTSJ3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn