

CSPP1 Rabbit pAb

货号: B15678

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	101kDa/141kDa/145kDa
实测分子量	145kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,K-562
细胞定位	Cytoplasm,centrosome,cytoskeleton,microtubule organizing center,spindle,spindle pole
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 962-1221 of human CSP P1 (NP_079066.5).
序列	QSPKGLDAATFQNVHDFNELKDRDSETRVDLKFMYLDPPRDHHTLEIQQQALLREQQKRLNRIKMQEGAKVDLDAIPSAK VREQRMPRDDTSDFLKNSLLESDSAFIGAYGETYPAIEDDVLPPPSQLPSARERRRNKWKGLDIDSSRPNVAPDGLSLKSI SSVNVDELVRVNEERMRLNEFHNPINTDDESSLVDPDDIMKHIGDDGSNSVATEPWLRPGTSETLKRFMAEQLNQEQ QQIPGKPGTFTWQGLSTAHG

靶点信息

研究背景	This gene encodes a centrosome and spindle pole associated protein. The encoded protein plays a role in cell-cycle progression and spindle organization, regulates cytokinesis, interacts with Nephrocystin 8 and i s required for cilia formation. Mutations in this gene result in primary cilia abnormalities and classical Joub ert syndrome. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
------	---

基因ID	79848
基因名	CSPP1
Swiss	Q1MSJ5
别名	CSPP1;CSPP;JBTS21

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn