

NRIP1 Rabbit pAb

货号: B15690

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	126kDa
实测分子量	150KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 939-1158 of human NRI P1 (NP_003480.2).
序列	LLLSENCVRDLSPHRSNSVADSKKKGHKNNVTNSKPEFSISLNLGLMYSSTQPSSCMDNRTFSYPGVVKTPVSPTFPEHLGCAGSRPESGLLNGCSMPSEKGPIKWVITDAEKNEYEKDSPRLTKTNPILYYMLQKGGNSVTSRETQDKDIWREASSAESVSQVTAKEELLPTAETKASFFNLRSPYNSHMGNNASRPHSANGEVYGLLGSVLTIKKESE

靶点信息

研究背景	Nuclear receptor interacting protein 1 (NRIP1) is a nuclear protein that specifically interacts with the hormone-dependent activation domain AF2 of nuclear receptors. Also known as RIP140, this protein modulates transcriptional activity of the estrogen receptor.
基因ID	8204

基因名	NRIP1
Swiss	P48552
别名	NRIP1;RIP140

产品验证



Western blot analysis of NRIP1 expressed in Mouse brain using NRIP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn