

RHOBTB1 Rabbit pAb

货号: B15693

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	79kDa
实测分子量	79kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa,Mouse kidney
细胞定位	cell cortex,cell projection,cytoplasmic vesicle,cytoskeleton,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 240-400 of human RHO BTB1 (NP_055651.1).
序列	PPVIKIPECPSMGTTNEAACLLDNPLCADVLFILQDQEHIFAHRIYLATSSSKFYDLFLMECEESPNGSEGACEKEKQSRDFQG RILSVDPEEEREEGPPRIPQADQWKSSNKSLEALGLEAEGAVPETQTLTGWSKGFIMHREMQRVNPISKRMGPMTVV

靶点信息

研究背景	The protein encoded by this gene belongs to the Rho family of the small GTPase superfamily. It contains a GTPase domain, a proline-rich region, a tandem of 2 BTB (broad complex, tramtrack, and bric-a-brac) domains, and a conserved C-terminal region. The protein plays a role in small GTPase-mediated signal transduction and the organization of the actin filament system. Alternate splicing results in multiple transcript variants.
基因ID	9886

基因名	RHOBTB1
Swiss	O94844
别名	RHOBTB1

产品验证



Western blot analysis of RHOBTB1 expressed in 293T,HeLa,Mouse kidney using RHOBTB1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn