

BAF60a Rabbit pAb

货号: **B15698**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	58kDa
实测分子量	58kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,293T,HepG2,BT-474,K-562,Mouse brain,Rat testis
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of mouse BAF60a (NP_114030.2).
序列	MAARAGFQSVAPSGGAGASGGAGVAAALGPGGTPGPPVRMGAPAGQGGLYRSPMPGAAYPRPGMLPGSRMTPQGPSM GPPGYGGNPSVRPGLAQSGMDQSRKRPAPQQIQVQQQAVQNRN

靶点信息

研究背景	The protein encoded by this gene is a member of the SWI/SNF family of proteins, whose members display helicase and ATPase activities and which are thought to regulate transcription of certain genes by alterin g the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent c hromatin remodeling complex SNF/SWI and contains a predicted leucine zipper motif typical of many tran scription factors.
------	---

基因ID	83797
基因名	
Swiss	Q96GM5
别名	CRACC1;Rsc8;SRG3;SWI3;SWI/SNF complex subunit SMARCC1;BAF60a;SMARCD1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn