

RRP7A Rabbit pAb

货号: B15777

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	32kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Rat heart
细胞定位	centrosome,cilium,cytoplasm,nucleolus,nucleoplasm,UTP-C complex
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human RRP7A (NP_056518.2).
序列	MVARRRKCAARDPEDRIPSPLGYAAIPIKFSEKQQASHYLYVRAHGVRQGTKSTWPQKRTLFLVNLNPPYCTEESLSRLLSTC GLVQSVELQEKPDLAESPKE SRSKFFHPKVPVGFQVAYVVFQKPSGVSAALALKGPLLVSTESHVPKSGIHKWISDYADSV PDPEALRVEVDTFMEAYDQKIAEEEEAKAKEEEGVPDEEGWVKVTRRGRRPVLP RTEAASLRVLERERRKRSRKELLNFYAW QHRESKMEHLAQLRKKFEEDKQRIELLRAQRKFRPY

靶点信息

研究背景	Enables RNA binding activity. Predicted to be involved in rRNA processing and ribosomal small subunit assembly. Predicted to act upstream of or within blastocyst formation. Predicted to be located in nucleoplasm. Predicted to be part of CURI complex and UTP-C complex. Implicated in primary autosomal recessive microcephaly.
------	--

基因ID	27341
基因名	RRP7A
Swiss	Q9Y3A4
别名	RRP7A;BK126B4.3;CGI-96;Rrp7

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn