

MEPE Rabbit pAb

货号: B15780

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	45kDa/58kDa/62kDa
实测分子量	58KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse liver,Mouse spinal cord,Rat brain,Rat liver
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 90-256 of human MEPE (NP_064588.1).
序列	FTNRQRLNKEYSISNKENTHNGLRMSIYPKSTGNKGFEDGDDAISKLHDQEEYGAALIRNNMQHIMGPVTAIKLLGEENKE NTPRNVLNIIIPASMNYAKAHSKDKKKPQRDSQAQKSPVKSSTHRIQHNIDYLKHLSKVKKIPSDFEGSGYTDLQERGDN DISPFS

靶点信息

研究背景	This gene encodes a secreted calcium-binding phosphoprotein that belongs to the small integrin-binding I igand, N-linked glycoprotein (SIBLING) family of proteins. Members of this family are components of the e xtracellular matrix of bone and dentin and regulate bone mineralization. Deficiency of a similar protein in mouse results in increased bone mass. Mice lacking this gene are resistant to aging-related trabecular bo ne loss. Alternative splicing results in multiple transcript variants.
------	---

基因ID	56955
基因名	MEPE
Swiss	Q9NQ76
别名	MEPE;OF45

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn