

PHLPP1 Rabbit pAb

货号: B15782

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	133kDa/184kDa
实测分子量	185KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,SH-SY5Y,Mouse testis,Mouse brain,Rat testis,Rat brain
细胞定位	Cell membrane,Cytoplasm,Membrane,Nucleus,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1160-1428 of human PH LPP1 (NP_919431.2).
序列	DQPSTGDASGAPAVWSHGYTEASGVKNKLCVAALSVNNFCDNREALYGVFDGDRNVEVPYLLQCTMSDILAEELQKTKN EEEYMVNTFIVMQRKLGTAGQKLGGAAVLCHIKHDPVDPGGSFTLTSANVGKQCQTVLCRNGKPLPLRSYIMSCEEELKRI KQHKAIITEDGKVNGVTESTRILGYTFLHPSVVPRPHVQSVLLTPQDEFFILGSKGLWDSLSVEEAVEAVRNPDALAAAKK LCTLAQSYGCHDSISAVVVQLSVTEDS

靶点信息

研究背景	This gene encodes a member of the serine/threonine phosphatase family. The encoded protein promotes apoptosis by dephosphorylating and inactivating the serine/threonine kinase Akt, and functions as a tumor suppressor in multiple types of cancer. Increased expression of this gene may also play a role in obesity and type 2 diabetes by interfering with Akt-mediated insulin signaling.
------	---

基因ID	23239
基因名	PHLPP1
Swiss	O60346
别名	PHLPP1;PHLPP;PLEKHE1;PPM3A;SCOP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn