

EYA1 Rabbit pAb

货号: **B15784**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:200 - 1:2000 IHC: 1:50 - 1:200
理论分子量	60kDa/61kDa/64kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,293T,Mouse testis,Mouse brain,Rat brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human EYA1 (NP_000494.2).
序列	NSLTNSSGFNSSQQDYPSYPSFGQGQYAQYYNSSPYAHYMTSSNTSPTTPSTNATYQLQEPPSGITSQAVTDPTAEYSTI HSPSTPIKDSDSDRLRRGSD

靶点信息

研究背景	This gene encodes a member of the eyes absent (EYA) family of proteins. The encoded protein may play a role in the developing kidney, branchial arches, eye, and ear. Mutations of this gene have been associated with branchiootorenal dysplasia syndrome, branchiootic syndrome, and sporadic cases of congenital cataracts and ocular anterior segment anomalies. A similar protein in mice can act as a transcriptional activator. Alternatively spliced transcript variants have been identified for this gene.
------	--

基因ID	2138
基因名	EYA1
Swiss	Q99502
别名	EYA1;BOP;BOR;BOS1;OFC1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn