

# ADAM28 Rabbit pAb

货号: B15788

产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                        |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  |                                                                                                        |
| 应用    | WB                                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                              |
| 理论分子量 | 61kDa/87kDa                                                                                            |
| 实测分子量 | 87kDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | Raji,HT-29,A-549,THP-1,Mouse liver,Rat liver                                                           |
| 细胞定位  | Cell membrane,Secreted,Single-pass type I membrane protein                                             |
| 纯化    | Affinity purification                                                                                  |

抗原信息


|      |                                                                                                                                                                                                                                                   |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 199-440 of human ADA M28 (NP_055080.2).                                                                                                                             |
| 序列   | VQEHEKYIEYYLVLDNGEFKRYNENQDEIRKRVFEMANYVNMLYKKLNTHVALVGMEIWTDKDKIKITPNASFTLENFSKW RGSVLSRRKRHDIAQLITATELAGTTVGLAFMSTMCSPPYSGVVQDHSNLLRVAGTMAHEMGMHNFHDDYSCKCP STICVMDKALSFYIPTDFSSCSRLSYDKFFEDKLSNCLFNAPLPTDIISTPICGNQLVEMGEDCDCGTSEECTNICCDAKTCK |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members o f this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, includin g fertilization, muscle development, and neurogenesis. The protein encoded by this gene is a lymphocyte -expressed ADAM protein. This gene is present in a gene cluster with other members of the ADAM family on chromosome 8. Alternative splicing results in multiple transcript variants. |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                          |
|-------|------------------------------------------|
| 基因ID  | 10863                                    |
| 基因名   | ADAM28                                   |
| Swiss | Q9UKQ2                                   |
| 别名    | ADAM28;ADAM 28;MDC-L;MDCL;eMDC II;eMDCII |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)