

# FBN2 Rabbit pAb

货号: B15803

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                          |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 157kDa/314kDa                                                                                        |
| 实测分子量 | 290kDa                                                                                               |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | LO2,Mouse brain                                                                                      |
| 细胞定位  | Secreted,extracellular matrix,extracellular space                                                    |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                                                                                                      |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 2733-2912 of human FB N2 (NP_001990.2).                                                                |
| 序列   | CVSGMGFNKGQYLSLDTEVDEENALSPEACYECKINGYSKKDSRQKRSIHEPDPTAVEQISLESVDMDSPVNMKFNLSHLGSKEHILELRPAIQPLNNHIRYVISQGNDDSVFRIHQRNGLSYLHTAKKKLMPGTYTLEITSIPLYKKKELKKLEESNEDDYLLGELGEALRMRLQIQLY |

## 靶点信息

|      |                                                                                                                                                                                                       |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene is a component of connective tissue microfibrils and may be involved in elastic fiber assembly. Mutations in this gene cause congenital contractural arachnodactyly. |
| 基因ID | 2201                                                                                                                                                                                                  |

|       |                   |
|-------|-------------------|
| 基因名   | FBN2              |
| Swiss | P35556            |
| 别名    | FBN2;CCA;DA9;EOMD |

产品验证



Western blot analysis of FBN2 expressed in LO2,Mouse brain using FBN2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)