

NFIX Rabbit pAb

货号: B15804

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	47kDa/48kDa/50kDa/54kDa/55kDa
实测分子量	45KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 180-260 of human NFIX (NP_002492.2).
序列	AYFVHTPESGQSDSSNQGDADIKPLPNGHLSFQDCFVTSGVWNVTELVRSQTPVATASGPNFSLADLESPSYNINQVT

靶点信息

研究背景	The protein encoded by this gene is a transcription factor that binds the palindromic sequence 5'-TTGGCN NNNNGCCAA-3 in viral and cellular promoters. The encoded protein can also stimulate adenovirus replication in vitro. Three transcript variants encoding different isoforms have been found for this gene.
基因ID	4784

基因名	NFIX
Swiss	Q14938
别名	NFIX;CTF;MRSHSS;NF-I/X;NF1-X;NF1A;SOTOS2

产品验证



Western blot analysis of NFIX expressed in A-549 using NFIX Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn