

EGFL7 Rabbit pAb

货号: B15806

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	29kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Mouse lung,Mouse kidney,Rat spleen,Rat kidney
细胞定位	Secreted,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 24-100 of human EGFL7 (NP_958854.1).
序列	YRPGRRVCAVRAHGDPVSESFVQRVYQPFLTTCDGHACSTYRTIYRTAYRRSPGLAPARPRYACCPGWKRTSGLPG

靶点信息

研究背景	This gene encodes a secreted endothelial cell protein that contains two epidermal growth factor-like domains. The encoded protein may play a role in regulating vasculogenesis. This protein may be involved in the growth and proliferation of tumor cells. Alternate splicing results in multiple transcript variants.
基因ID	51162
基因名	EGFL7

Swiss	Q9UHF1
别名	EGFL7;NEU1;VE-STATIN;ZNEU1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn