

ZC3HAV1 Rabbit pAb

货号: B15824

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	18kDa/42kDa/67kDa/77kDa/101kDa
实测分子量	101kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,LO2,OVCAR3,Mouse lung,Mouse thymus
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 150-400 of human ZC3 HAV1 (NP_078901.3).
序列	CKSYKGEGRQQICNQPPCSRLHICDHFTRGNCRFPNCLRSHNLMDRKVLAIMREHGLNPDVVQNIQDICNSKHMQKNP PGPRAPSSHRRNMAYRARSKSRDRFFQGSQEFLASASASAERSCTPSPDQISHRASLEDAPVDDLTRKFTYLGSQDRAR PPSGSSKATDLGGTSQAGTSQRFLENGSQLDLLHGPNPGSTYLASNSTSAPNWKSLTSWTNDQGARRKTVFSPTLPAARS SLGSLQTPEAVTTR

靶点信息

研究背景	This gene encodes a CCCH-type zinc finger protein that is thought to prevent infection by retroviruses. Studies of the rat homolog indicate that the protein may primarily function to inhibit viral gene expression and induce an innate immunity to viral infection. Alternative splicing occurs at this locus and two variants, each encoding distinct isoforms, are described.
------	--

基因ID	56829
基因名	ZC3HAV1
Swiss	Q7Z2W4
别名	ZC3HAV1;ARTD13;FLB6421;PARP13;ZAP;ZC3H2;ZC3HDC2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn