

NUP160 Rabbit pAb

货号: B15837

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:20 - 1:50
理论分子量	24kDa/55kDa/162kDa
实测分子量	160kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-431,Jurkat,HepG2
细胞定位	Nucleus,nuclear pore complex
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 40-150 of human NUP160 (NP_056046.1).
序列	ALERSFVELSGAERERPRHFREFTVCSIGTANAVAGAVKYESAGGFYYVESGKLFSVTRNRFIHWKTSGDTLELMEESLDI NLLNNAIRLKFQNCVLPGGVYVSETQNR

靶点信息

研究背景	NUP160 is 1 of up to 60 proteins that make up the 120-MD nuclear pore complex, which mediates nucleoplasmic transport.[supplied by OMIM, Apr 2004]
基因ID	23279

基因名	NUP160
Swiss	Q12769
别名	NUP160

产品验证



Western blot analysis of NUP160 expressed in HeLa,A-431,Jurkat using NUP160 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn