

# TCF7L1 Rabbit pAb

货号: B15858

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000                                                                            |
| 理论分子量 | 62kDa                                                                                                |
| 实测分子量 | 75KDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | HeLa,HCT116,A-549,A-431,Mouse lung,Mouse brain,Mouse skin,Rat lung,Rat skin                          |
| 细胞定位  | Nucleus                                                                                              |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                       |
|------|-----------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 260-340 of human TCF7 L1 (NP_112573.1). |
| 序列   | YSLPPGGFRHPYPALAMNASMSSLVSSRFSPHMVAPAHPLPTSGIPHPAIVSPIVKQEPAPPSLSPAVSVKSPVTVKKEE                                      |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the T cell factor/lymphoid enhancer factor family of transcription factors. These transcription factors are activated by beta catenin, mediate the Wnt signaling pathway and are ant agonized by the transforming growth factor beta signaling pathway. The encoded protein contains a high mobility group-box DNA binding domain and participates in the regulation of cell cycle genes and cellular senescence. |
| 基因ID | 83439                                                                                                                                                                                                                                                                                                                                                                                                                                           |

|       |                   |
|-------|-------------------|
| 基因名   | TCF7L1            |
| Swiss | Q9HCS4            |
| 别名    | TCF7L1;TCF-3;TCF3 |

产品验证



Western blot analysis of TCF7L1 expressed in Mouse lung,Rat lung,HeLa using TCF7L1 Rabbit pAb at 1: 1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)